

## EndoView CH Series Single use Flexible Choledochoscope



### Technical Specifications

Model		CH-M32	CH-M40	CH-M50	CH-M52
Optical system	Field view	120°			
	Depth of field	3-50 mm			
Insertion tube	Distal end outer diameter	3.2mm	4.4mm	5.0mm	5.2mm
	Working length	350 mm			
Working Channel		1.2mm	> 1.1mm	2.8mm	> 2.1mm
Angulation	Bending	210 ° U/D			
Total Length		610mm			

Processor	
Function	Photo
	Video
	White Balance
	Freeze
Storage	Minimal 16G
Connection	DVI, HDMI
	VGA
	SDI, AV

### Shenzhen HugeMed Medical Technical Development Co.,Ltd

Add: 416-1,516-1, Building 2, No. 1, Mawu Road, Baoan Community, Yuanshan Street, Longgang District, Shenzhen 518115, China  
Tel: +86-755-2227 5866  
+86-755-8226 7820  
Fax: +86-755-2227 5833  
Email: sale@hugemed.net

[www.hugemed.net](http://www.hugemed.net)

## AN ENDOSCOPE SOLUTION COMPANY

EndoView CH Series  
Single use Flexible  
Choledochoscope

Always Available  
Always Brand New  
No Costly Repair  
No Contamination

# EndoView CH Series Single use Flexible Choledochoscope

- All in one system, plug and play
- Always available, avoiding procedure delay



210° / 210°

high deflection, full  
navigate

Compatible with conventional  
accessories



Portable and light

Customized Two buttons: WB /  
Photo / Video / Freeze

Standard **1 Million** Pixel  
High Quality Image