

EndoView BR Series Single use Flexible Bronchoscope



Technical Specifications

Model		BR-M22	BR-M32	BR-M40	BR-M50	BR-M52	BR-M58	BR-M58B
Optical system	Field view	120°						
	Depth of field	3-50mm						
Insertion tube	Distal end outer diameter	2.2mm	3.2mm	4.4mm	5.0mm	5.2mm	5.8mm	5.8mm
	Working length	600mm						
Working Channel		/	1.2mm	2.0mm	2.8mm	2.2mm	2.8mm	3.3mm
Angulation	Bending	210° U/D						
Total Length		860mm						

Processor	
Function	Photo
	Video
	White Balance
	Freeze
Storage	Minimal 16G
Connection	HDMI
	SDI, AV
Display	
Screen	LCD
	Size: 8 inch
Function	Photo
	Video
	White Balance
	Freeze
Battery	Type: Lithium
	Maximum 4 hours
Connection	HDMI

AN ENDOSCOPE SOLUTION COMPANY

EndoView BR Series
Single use Flexible Bronchoscope

- Always Available
- Always Brand New
- No Costly Repair
- No Contamination




Shenzhen HugeMed Medical Technical Development Co.,Ltd

Add: 416-1,516-1, Building 2, No. 1, Mawu Road, Baoan Community, Yuanshan Street, Longgang District, Shenzhen 518115, China
Tel: +86-755-2227 5833
Fax: +86-755-2227 5833
Email: info@hugemed.net

www.hugemed.net


EndoView BR Series Single use Flexible Bronchoscope

- Plug in and ready to go
- Multiple size to satisfy various clinical settings
- Navigate with ease



Standard **1 Million** Pixels,
High Quality Image

- Slim: OD 3.2/ ID 1.2 mm
- Regular: OD 4.4/ ID 2.0 mm
- Regular: OD 5.0/ ID 2.8 mm
- Large: OD 5.8/ ID 3.4 mm



210° / 210° high deflection,
Ensure full navigation



Customized Two buttons:
WB / Photo / Video / Freeze



USB and HDMI connection,
easy data transfer



Compatible with conventional
accessories